

**AVAILABLE FOR TRIAL
AT OUR PARTNER
LOCATIONS**



"I truly realized the relief Skeed provides for my hands, wrists, and arms. You just activate the lever with your fingers and it brakes! It's a daily advantage. Even indoors or at low speeds, nothing compares!"

Clément
Skeed user

**TO TEST SKEED OR BEGIN
A PARTNERSHIP:**

✉ hello@eppur.eu 📞 +33.6.84.80.93.33

WWW.EPPUR.EU



**SCAN HERE TO
CONTACT US:**



Do not dispose of in public areas.

EPPUR

-SKEED

GROUND BRAKING TECHNOLOGY



MADE IN FRANCE

This medical device is a regulated health product which, under this regulation, bears the CE marking.

EPPUR, a simplified joint-stock company with a share capital of €28,901, registered with the Registry of Commerce and Companies of Lille Métropole under number B 881 789 408.

STEER AND BRAKE

WITH YOUR FINGERTIPS

Skeed is a pair of wheels compatible with all manual wheelchairs, featuring a **unique braking system built into the hub**.

A fingertip-accessible brake lever is connected to the wheel's center to deliver **smooth, progressive braking**.



LEVER GRIP
TETRA OPTION



LEVER GRIP
PARA OPTION



PROGRESSIVE BRAKE

The more you pull, the more you brake !
Activate the brake lever backwards to slow down and steer your wheelchair with unprecedented precision.



PARKING BRAKE

Push the lever forward to immobilize your wheelchair. A powerful brake ensures total stability, regardless of your tire pressure.

Ask for our order sheet model to check all the products options and details.

LESS EFFORT, MORE SAFETY

SKEED SAVES YOUR HANDS AND ARMS



100% WHEELIE CONTROL PRESERVED

Thanks to its fixed handrim, Skeed guarantees you perfect wheelie control.

Overcome obstacles and maneuver with agility.



5X LESS EFFORT AND ZERO FRICTION

With Skeed, no more burns and risk of injury ! The lever eliminates all friction on the hands and divides the braking effort by 5.

Your arms and hands are thus preserved !



COMPATIBLE WITH MOTORIZATION

Skeed wheels are especially useful for mastering fifth rear propulsion wheels safely, improving steering control without straining the shoulders on every turn or slope.



ULTRA-LIGHT FULLY EQUIPPED WHEEL

Only 2.15 Kg* to facilitate transfers and reduce the overall weight of your wheelchair !

*Standard configuration
Naked rim, hub & lever : 1.5 Kg